



# HESSBURG-BARRON VACUUM TREPHINE

with uninterrupted 360° suction chamber



**ORIGINAL DESIGN WITH  
OVER 20 YEARS  
OF PROVEN RESULTS**

- crosshairs assist with centration
- self contained with its own suction source
- trephine of choice worldwide
- pre-sterile single use
- guaranteed sharp




**Maintains a Perpendicular Relationship to the Cornea**

# Hessburg-Barron Vacuum Trepine

The Hessburg-Barron trephine provides complete or partial keratotomy for keratoplasty procedures.

Unlike conventional units this vacuum trephine attaches to the cornea during the trephining process to hold the cornea in its natural state allowing for a perfect perpendicular cut for the recipient well.



The trephine incorporates cross hairs to aid in proper centration over the optical axis.

- 21-8260 Hessburg-Barron Vacuum Trepine, 6.0mm
- 21-8265 Hessburg-Barron Vacuum Trepine, 6.5mm
- 21-8270 Hessburg-Barron Vacuum Trepine, 7.0mm
- 21-8272 Hessburg-Barron Vacuum Trepine, 7.25mm
- 21-8275 Hessburg-Barron Vacuum Trepine, 7.5mm

- 21-8277 Hessburg-Barron Vacuum Trepine, 7.75mm
- 21-8280 Hessburg-Barron Vacuum Trepine, 8.0mm
- 21-8282 Hessburg-Barron Vacuum Trepine, 8.25mm
- 21-8285 Hessburg-Barron Vacuum Trepine, 8.5mm
- 21-8287 Hessburg-Barron Vacuum Trepine, 8.75mm
- 21-8290 Hessburg-Barron Vacuum Trepine, 9.0mm

## Barron Marking Corneal Punch



The Barron disposable punch system is designed to make an exact cut of the donor tissue for corneal transplants. The ultra-sharp cutting blade is guided by four posts to assure proper alignment. After the tissue is properly placed in the punch, an ink pen is used to mark the donor button in four symmetric locations. These four marks will aid in the placing of primary sutures.

The system is complete with punch, blade, block and marking pen, sterile packed and ready for use. Because this system is totally disposable, you are assured a clean straight punch each time.

- 21-8360 Barron Marking Corneal Punch System, 6.00mm
- 21-8365 Barron Marking Corneal Punch System, 6.50mm
- 21-8367 Barron Marking Corneal Punch System, 6.75mm
- 21-8370 Barron Marking Corneal Punch System, 7.00mm
- 21-8372 Barron Marking Corneal Punch System, 7.25mm
- 21-8375 Barron Marking Corneal Punch System, 7.50mm
- 21-8377 Barron Marking Corneal Punch System, 7.75mm
- 21-8380 Barron Marking Corneal Punch System, 8.00mm
- 21-8382 Barron Marking Corneal Punch System, 8.25mm
- 21-8385 Barron Marking Corneal Punch System, 8.50mm
- 21-8387 Barron Marking Corneal Punch System, 8.75mm
- 21-8390 Barron Marking Corneal Punch System, 9.00mm
- 21-8395 Barron Marking Corneal Punch System, 9.50mm



DISTRIBUIDO POR



Concilio de Trento, 37-39  
08018 Barcelona  
Tel. 93 303 65 65  
bloss@blossgroup.com  
www.blossgroup.com