





**Lenstec
Capsular Rings...
The Blue
Makes It
a Stand Out!**

**Acrylic Foldable,
Silicone Foldable and
PMMA IOLs**

**LENS
TEC** 

Capsular Rings



Model No.	Size		Color
	Open	Compressed	
LR-1300	12.30	10.00	Clear
LR-1400	14.30	12.00	Clear
LR-1300B	12.30	10.00	Blue
LR-1400B	14.30	12.00	Blue

Over ten years ago the engineers at Lenstec developed the "visible" Capsular Tension Ring so that surgeons can focus on the implantation procedure rather than looking for the ring.

Type LR-1300B and LR-1400B are made of blue PMMA material. They offer high flexibility, high memory control and are easy to insert manually or with any available injector.

Both designs are applicable when the zonules are missing or damaged, in cases of lens luxation, during phaco

complications, zonular breakage during sudden unfolding of a foldable lens, for repositioning of a lens or if decentration of same. Other uses apply for myopia, pseudoexfoliation, eyes undergoing vitrectomy and trauma.

The implantation of each design offers various advantages: force distribution acting on the cataract and capsule circumferentially along the entire device, less risk of capsular shrinking, less risk of capsular fibrosis, and safer IOL centration.



World Headquarters

Lenstec, Inc.
1765 Commerce Ave. N.
St. Petersburg, FL 33716 USA
Tel: 727-571-2272 • Fax: 727-571-1792
Web: www.lenstec.com
Email: lenstec@lenstec.com

Lenstec Barbados

Lenstec Barbados Inc.
Airport Commercial Centre
Pilgrim Road
Christ Church, Barbados BB17092
Tel: 246-420-6795 • Fax: 246-420-6797
Email: lenstecbarbados@lenstec.com

DISTRIBUIDO POR

Concilio de Trento, 37-39
08018 Barcelona
Tel. 93 303 65 65
bloss@blossgroup.com
www.blossgroup.com